




P(LA/CL)

Monofilament Bioabsorbable Suture Material

P(LA/CL) monofilament bioabsorbable suture material has a good reputation among surgeons for being easy to handle. The suture loses 50% of its initial tensile strength in about 8 weeks according to In-vitro degradation tests. Our **P(LA/CL)** is currently available in many countries overseas. The sterilized finished product of our **P(LA/CL)** cannot only be used for human surgical procedures, but also for veterinary procedures as well.

Composition

P(LA/CL): Poly (L-lactide-co- ε - caprolactone)



Characteristics

- Low incidence of infection
- Superior smoothness
- Superior tensile strength

Product Variation

Color: Violet and Beige

Size: USP2 to USP8-0 are available



Retention of Tensile strength

